

Scrub Stick Dentifrice

Formulated with Omyacare® and Omyadent®, this waterless toothpaste has everything you need to clean and remineralize teeth in a revolutionary new format.



Scrub Stick Dentifrice



Ingredients	INCI Nomenclature	Suppliers	% w / w
A KahlWax 2039L Candelilla	Euphorbia Cerifera (Candelilla) Wax	KahlWax	10.20
Resplanta® AVOCADO MB	Avocado Oil Glycereth-8 Esters, Water	Res Pharma Industriale	7.00
Coconut Oil	Cocos Nucifera (Coconut) Oil	---	20.00
Jobaba Oil	Simmondsia Chinensis (Jojoba) Seed Oil	---	15.00
Omyacare® S 70-AZ	Calcium Carbonate	Omya	40.00
Omyadent® 100-OG	Hydroxyapatite, Calcium Carbonate	Omya	5.00
Erylite® Stevia 400	Erythritol, Stevia Rebaudiana Extract	Jungbunzlauer	0.80
Peppermint Oil	Mentha Piperita Oil	---	2.00

100.00

Procedure

1. Combine all ingredients. Heat to 80°C and mix until homogeneous.
2. Pour mixture into molds and cool to room temperature.

How To Use

Rub toothbrush bristles on dentifrice stick and brush teeth for two minutes.



Omya Specialty Materials Inc.
5150 E. Pacific Coast Hwy, Suite 600
Long Beach, CA 90804
(562) 961-3333 | PersonalCareUS@omya.com

Innovation Comes Naturally™

Scan To Learn More

